

View Protocols SPA

Edit Copy Delete Steps Browse Db Collapse Labels

- ▼ Picking
 - reion - auto-picking
 - xmipp3 - particle boxsize
 - xmipp3 - manual-picking (step 1)
 - xmipp3 - auto-picking (step 2)
 - xmipp3 - deep consensus picking
 - xmipp3 - deep micrograph cleaner
 - xmipp3 - assign tiltpairs
 - xmipp3 - picking consensus
 - xmipp3 - pick noise
 - sphire - cryolo picking
 - sphire - cryolo training
 - eman2 - boxer
 - eman2 - boxer auto
 - eman2 - sparx gaussian picker
 - bsoft - particle picking

- ▼ Extract
 - reion - particles extraction
 - reion - movie particles extraction
 - reion - export particles
 - reion - sort particles
 - xmipp3 - extract particles

- ▶ more
- ▶ Preprocess
- ▶ Filter
- ▶ Mask

- ▼ 2D
 - Align
 - Classify

- ▼ 3D
 - Initial volume
 - Preprocess
 - Classify
 - Refine
 - Postprocess

Protocol: xmipp3 - extract particles
finished
Cite
Help

Run

Run name xmipp3 - extract particles smal **Comment**

Run mode ☒ Continue ☐ Restart ? **Host** localhost

Parallel ☒ Threads ☐ MPI 4 ? **Use queue?** ☐ Yes ☒ No ?

Wait for ?

Expert Level ☒ Normal ☐ Advanced

Input Preprocess

Input

Input coordinates xmipp3 - deep micrograph cleaner.outputCoordinates_Manual_096 ?

Micrographs source ☒ same as picking ☐ other ?

CTF estimation xmipp3 - ctf consensus.outputCTF. ?

Downsampling factor 4.0 ?

Particle box size (px) 74 ?

Fill pixels outside borders ☐ Yes ☒ No ?

Close Save Execute


Input

- inputCoordinates (from xmipp3 - deep micrograph cleaner -> outputCoordinates Manuz SetOfCoordinates (6431 items, 300 x 300)
- ctfRelations (from xmipp3 - ctf consensus -> outputCTF) SetOfCTF (47 items)

Output

- xmipp3 - extract particles small -> outputParticles SetOfParticles (6368 items, 74 x 74, 1.98 Å/px)

View: Tree Refresh



Analyze Results