

SCIPION

v3.0.6 - Eugenius

View Protocols SPA

- Protocols SPA
 - Imports
 - import movies
 - import micrographs
 - import particles
 - import volumes
 - more
 - Movies
 - Micrographs
 - CTF estimation
 - Preprocess
 - Particles
 - Picking
 - Extract
 - Preprocess
 - Filter
 - Mask
 - 2D
 - Align
 - Classifu
 - 3D
 - Initial volume
 - Preprocess
 - Classifu
 - Refine
 - Postprocess
 - xmipp3 - localdeblur sharpening**
 - xmipp3 - deepEMhancer
 - relion - post-processing
 - Analysis
 - Reconstruct
 - Tools
 - Exports

Edit Copy Delete Browse

xmipp3 - highres angles&shifts 2, outputVolume finished

xmipp3 - create 3d mask finished

noD... xmipp3 - local deepRes finished

xmipp3 - crop/resize volumes resolution, outputVolume finished

xmipp3 - localdeblur sharpening finished

Summary Methods Output Log


Input

- inputVolume (from xmipp3 - highres angles&shifts 2 -> outputVolume) Volume (250 x 250 x 250, 0.74 Å/px)
- resolutionVolume (from xmipp3 - crop/resize volumes resolution -> Volume (250 x 250 x 250, 0.74 Å/px)

Output

- SUMMARY
- LocalDeblur Map

Protocol Run: XmippProtLocSharp



Protocol: xmipp3 - localdeblur sharpening

finished Cite Help

Run

Run name xmipp3 - localdeblur sharpening Comment

Run mode ☒ Continue ☐ Restart ? Host localhost

Parallel Threads 4 ? Use queue? ☐ Yes ☒ No ?

Wait for

Expert Level ☒ Normal ☐ Advanced

Input

Input Map xmipp3 - highres angles&shifts 2, outputVolume ?

Resolution Map xmipp3 - crop/resize volumes resolution, outputVolume ?

Close Save Execute